

Vapour Pressure Osmometry & Light scattering Technique

Vapor-Phase Osmometry



Based on vapour pressure lowering: Colligative property

 $\overline{M_n}$

P₁ is vapour Pressure of solvent:

If solute is non volatile:

After rearranging:

ΔP: lowering in vapour Pr

$$P = P_1 + P_2 = x_1 P_1^0 + x_2 P_2^0$$

$$P = x_1 P_1^0 = (1 - x_2) P_1^0$$

$$x_2 = 1 - \frac{P}{P_1^0} = \frac{P_1^0 - P}{P_1^0} = \frac{-\Delta P}{P_1^0}$$

$$\frac{\Delta P}{P_1^0} = -x_2 = \frac{W_2/M_2}{(W_1/M_1) + (W_2/M_2)}$$

$$\lim \left(\frac{\Delta P}{P_1^0}\right)_{C_1 \to 0} = -\frac{(W_2/M_2)}{(W_1/M_1)} = -\frac{W_2M_1}{W_1M_2} = -c_2\frac{V_1^o}{M_2}$$

Vapor-Phase Osmometry contd...



Writing in virial form:

$$\frac{\Delta P}{P_1^0} = -V_1^0 c_2 \left[\frac{1}{M} + B c_2 + C c_2^2 + \cdots \right]$$

We can relate ΔP with ΔT using Clausius Clapeyron equation:

$$\Delta T = \frac{\Delta P}{P_1^0} \cdot \frac{RT^2}{\Delta H_v}$$
Latent heat of vaporization of solvent

$$\frac{\Delta T}{c_2} = -\frac{RT^2}{\Delta H_v} V_1^0 \left[\frac{1}{M} + Bc_2 + Cc_2^2 + \cdots \right]$$

 M_n of solute can now estimate by measuring $\frac{\Delta T}{c_2}$ while extrapolating it to $c_2=0$

Practical Aspects




- NO need of membrane in VPO.
- k_s is an instrument constant
- Depend on T, solvent, thermistor pair

$$\frac{\Delta T}{c_2} = k_s \left[\frac{1}{\overline{M}_n} + Bc_2 + Cc_2^2 + \cdots \right]$$

\[
 \Delta T \] can be measured very accurately as a fn of the bridge imbalance output voltage

$$\frac{\Delta V}{c} = \frac{K}{M} + KBc$$



Schematic diagram of a vapour phase osmometer

$$K = M(\Delta V/c)_{c\to 0}$$


$$\longrightarrow \text{Calibration constant}$$

Illustration



VPO data for PS of known MW & HTPB in Toluene (70°C). MW of HTPB...???

Polymer	Concentration, c (g/L)	Bridge output, ΔV (μV)	
Standard	6	107	ΔV K
polystyrene	9	164	$\frac{\Delta V}{} = \frac{K}{R} + KBc$
of $\overline{M}_n = 1800$	12	224	$\frac{-1}{c} = \frac{-1}{M} + KBc$
	15	28 7	
HTPB of unknown	6	85	
molecular weight	9 ,	129	$K = M(\Delta V/c)_{c\to 0}$
· ·	12	176	
	15	225	



Light-Scattering Method



- Important technique for determination of Mw.
- Provide information about size and shape of polymer molecules in solution
- Also characterize the interaction between solvent and polymer molecules
- Large difference in intensity of incident beam and light scattered beam from polymer solution.


Rayleigh Ratio (R_{θ}):

Equation for gaseous molecules:

$$\frac{i_{\theta} r^2}{I_0(1 + \cos^2 \theta)} = R_{\theta} = \frac{2\pi^2}{\lambda^4 N_{\text{Av}}} \frac{(\widetilde{n} - 1)^2 M}{c}$$

c: conc. of gas (mass/vol)

Can not applicable to the liquids due to strong intermolecular forces are present in Liquids



Arrangement of apparatus to measure light scattered from solution

Light-Scattering Method contd...



Einstein developed approach for liquids:

Based on local fluctuations in the density, so in refractive index and hence in scattering of incident light.

$$R_{\theta} = \frac{i_{\theta} r^2}{I_0 (1 + \cos^2 \theta)} = \frac{2\pi^2}{\lambda^4 N_{\text{Av}}} \frac{RT}{\beta} \left(\tilde{n} \frac{d\tilde{n}}{dp} \right)^2$$

$$R_{\theta}(\text{solution}) - R_{\theta}(\text{solvent})$$

Presence of solute in solution causes additional scattering of light (Debye model incorporate local fluctuation due to solute concentration:

Random thermal motion are opposed by the osmotic pressure of the solution:

Excess scattering

from liquid due to the solute
$$R'_{\theta} = \frac{i'_{\theta} r^2}{I_0(1 + \cos^2 \theta)} = \frac{2\pi^2}{\lambda^4 N_{\text{Av}}} \left(\tilde{n}_0 \frac{d\tilde{n}}{dc}\right)^2 \frac{RTc}{(\partial \Pi/\partial c)_T}$$

 $R_{\theta}(\text{solution}) - R_{\theta}(\text{solvent})$ Specific RI increment with c


Light-Scattering Method contd...



❖ Osmotic pressure for monodisperse solute

$$\left(\frac{\partial\Pi}{\partial c}\right)_T = RT\left(\frac{1}{M} + 2A_2c + 3A_3c^2 + \cdots\right)$$

❖Can be estimate by galvanometer readings (I) directly



Relation b/w Turbidity and Rayleigh Ratio


- **❖Turbidity:** is measure of decrease of the incident beam Intensity (I₀) scattered I_s per unit length, after passing through of solution.
 - Can be related as bear's law:

$$I/I_0 = e^{-\tau \ell}$$


$$e^{-\tau \ell} = \frac{I_0 - I_s}{I_0} = \frac{I_0 - I_s' \ell}{I_0} = 1 - \frac{I_s' \ell}{I_0}$$
series:
$$\tau = I_s' / I_0$$

$$e^{-\tau \ell} = 1 - \tau \ell + \frac{1}{2} (\tau \ell)^2 - \frac{1}{6} (\tau \ell)^3 + \cdots \simeq 1 - \tau \ell$$

After expanding this series:



$$-\tau' = 1 - \tau' + \frac{1}{2}(\tau')^2 - \frac{1}{6}(\tau')^3 + \cdots \simeq 1 - \tau'$$



The variation of effective scattering volume (shaded) with angle:

(a) $\theta < 90^{\circ}$; (b) $\theta = 90^{\circ}$. A: incident beam; B: polymer solution in solvent; C: scattering cell; D: light trap; E: scattered light beam; F: movable photocell connected to galvanometer.

Polymer Science & Technology (CL-623)

Relation b/w Turbidity and Rayleigh Ratio

Total light intensity scattered per unit length through all angles polar coordinates

$$I_s' = \int_0^{\pi} \int_0^{2\pi} r^2 i_{\theta}' \sin\theta d\theta d\phi$$

$$R'_{\theta} = \frac{i'_{\theta}r^2}{I_0(1+\cos^2\theta)}$$

$$\frac{I_s'}{I_0} = \tau = \int_0^{\pi} \int_0^{2\pi} R_{\theta}'(1 + \cos^2\theta) \sin\theta d\theta d\phi$$

$$R'_{\theta} = \frac{2\pi^2}{\lambda^4 N_{\text{Av}}} \left(\tilde{n}_0 \frac{d\tilde{n}}{dc} \right)^2 \frac{RTc}{(\partial \Pi/\partial c)_T}$$

$$\tau = \frac{2\pi^2}{\lambda^4 N_{\text{Av}}} \left(\tilde{n}_o \frac{d\tilde{n}}{dc} \right)^2 \frac{RTc}{(\partial \Pi/\partial c)_T} \int_0^{\pi} (1 + \cos^2 \theta) \sin \theta d\theta \int_0^{2\pi} d\phi$$

$$\tau = \frac{I_s'}{I_0} = \left(\frac{32\pi^3}{3\lambda^4 N_{\text{AV}}}\right) \frac{\left(\tilde{n}_o \frac{d\tilde{n}}{dc}\right)^2 RTc}{(\partial \Pi/\partial c)_T}$$

$$\tau = \frac{16\pi}{3} R_\theta' = \left(\frac{16\pi}{3}\right) \frac{i_\theta' r^2}{I_0(1 + \cos^2 \theta)}$$

$$16\pi/3$$

$$\tau = \frac{16\pi}{3}R'_{\theta} = \left(\frac{16\pi}{3}\right)\frac{i'_{\theta}r^2}{I_0(1+\cos^2\theta)}$$

Turbidity and Molecular Weight of Polymer



Light scattering of solution may now be related to the solute molecular weight

$$\frac{\Pi}{c} = RT \left(\frac{1}{M} + A_2 c + A_3 c^2 + \cdots \right) \qquad \left(\frac{\partial \Pi}{\partial c} \right)_T = RT \left(\frac{1}{M} + 2A_2 c + 3A_3 c^2 + \cdots \right)$$

$$\tau = \frac{I_s'}{I_0} = \left(\frac{32\pi^3}{3\lambda^4 N_{\text{Av}}} \right) \frac{\left(\tilde{n}_0 \frac{d\tilde{n}}{dc} \right)^2 RT c}{(\partial \Pi/\partial c)_T} \qquad \Rightarrow \tau = \left(\frac{32\pi^3}{3\lambda^4 N_{\text{Av}}} \right) c \left(\tilde{n}_0 \frac{d\tilde{n}}{dc} \right)^2 / \left(\frac{1}{M} + 2A_2 c + 3A_3 c^2 + \cdots \right)$$

$$H = \left(\frac{32\pi^3}{3\lambda^4 N_{\rm Av}}\right) \left(\tilde{n}_0 \frac{d\tilde{n}}{dc}\right)^2$$


$$\frac{Hc}{\tau} = \frac{1}{M} + 2A_2c + 3A_3c^2 + \cdots$$

$$K = \frac{3H}{16\pi} = \left(\frac{2\pi^2}{\lambda^4 N_{AV}}\right) \left(\tilde{n}_0 \frac{d\tilde{n}}{dc}\right)^2 \qquad \frac{Kc}{R'_{\theta}} = \frac{3Hc}{16\pi R'_{\theta}} = \frac{Hc}{\tau} = \frac{1}{M} + 2A_2^2c + 3A_3c^2 + \cdots$$




Light scattering determination of the molecular weight of polystyrene

$$\frac{Hc}{\tau} = \frac{1}{M} + 2A_2c + 3A_3c^2 + \cdots$$



$$\frac{Kc}{R'_{\theta}} = \frac{3Hc}{16\pi R'_{\theta}} = \frac{Hc}{\tau} = \frac{1}{M} + 2A_{2}c + 3A_{3}c^{2} + \cdots$$

Plot of K_c/R'₉₀ versus c in the absence of interference effect



Example:



- Typical light scattering data for solutions of polystyrene in benzene (n_0 =1.5010) with λ =4358Å at 20°C
- ❖ Determine the molecular weight of the polymer?

$$H = \left(\frac{32\pi^3}{3\lambda^4 N_{\text{Av}}}\right) \left(\tilde{n}_0 \frac{d\tilde{n}}{dc}\right)^2$$

$$H = \left(\frac{32\pi^3}{3\lambda^4 N_{\text{Av}}}\right) \tilde{n}_0^2 \left(\frac{\tilde{n} - \tilde{n}_0}{c}\right)^2$$


$$\left(\frac{Hc}{T}\right)_{c=0}$$

Conc. (g/100 cm ³)	$(\tilde{n}-\tilde{n}_o)\times 10^2$	$\tau \times 10^{2} (\text{cm}^{-1})$
0.175	0.021	0.093
0.385	0.043	0.144
0.594	0.066	0.182
0.730	0.082	0.198
1.000	0.110	0.226

$$\frac{Hc}{\tau} = \frac{1}{M} + 2A_2c + 3A_3c^2 + \cdots$$

Dissymmetry of Scattering







$$z_{\beta} = \frac{i'_{90}^{\circ} - \beta}{i'_{90}^{\circ} + \beta} = \frac{i'_{\theta}}{i'_{\pi - \theta}} = \frac{R'_{\theta}}{R'_{\pi - \theta}}$$

$$P(\theta) = \frac{(R'_{\theta})_{\text{observed}}}{(R'_{\theta})_{\text{no interference}}} = \frac{(i'_{\theta})_{\text{observed}}}{(i'_{\theta})_{\text{no interference}}}$$



$$\frac{Kc}{(R'_{\theta})_{\text{observed}}} = \frac{Kc}{P(\theta)(R'_{\theta})_{\text{no interference}}}$$
$$= \frac{1}{P_{\theta}} \left[\frac{1}{M} + 2A_{2}c + 3A_{3}c^{2} + \cdots \right]$$

$$\frac{Kc}{R'_{\theta}} = \frac{Hc}{\tau} = \frac{1}{M.P(\theta)} + 2A_2c$$

$$\left(\frac{Kc}{R'_{\theta}}\right)_{c\to 0} = \left(\frac{Hc}{\tau}\right)_{c\to 0, \ \theta\to 0} = \frac{1}{M}$$

Scattering Factor



$$P(\theta) = \frac{2}{u^2} \left[e^{-u} - (1-u) \right]$$

$$u = \left(\frac{2}{3}\right) \left(\frac{\langle r^2 \rangle}{{\lambda'}^2}\right) \left(2\pi \sin \frac{\theta}{2}\right)^2$$


$$P(\theta) = 1 - u/3$$

$$1/P(\theta) = 1 + \left(\frac{8\pi^2}{9{\lambda'}^2}\right) \langle r^2 \rangle \sin^2 \frac{\theta}{2}$$



$$\frac{Kc}{R'_{\theta}} = \frac{Hc}{\tau} = \frac{1}{M} + \left(\frac{1}{M}\right) \left(\frac{8\pi^2}{9\lambda'^2}\right) \langle r^2 \rangle \sin^2(\theta/2) + 2A_2c$$

Zimm Plots



Zimm plot for the polystyrene sample